

Oxygen Monitor



FEATURES:

- Fast Response.
- Accurate & Reliable.
- Auto calibration for 21% & 100% oxygen.
- Operated on Mains supply.
- User selectable audiovisual alarm limits.
- Both upper & lower alarm limits selectable by only two key.
- Large backlight LCD.
- Size in cm 17 x 6 x 13.
- Battery backup (optional).

Technical Specification:

Ambient operating temperature Range:	10°C to 50°C.
Operating Altitude:	-557 ft. to +16,000 ft / -170m to +4877m.
Monitoring range:	18% to 100%.
Sensor used:	Galvanic sensor.
Transport/Storage Temperature:	-13°F to +131°F (-25°C to +55°C).
Accuracy:	+/- 1 %.
Ambient working Humidity:	15% to 95% non-condensing.
Dimensions:	(17cm x 6cm x 13cm).
Weight:	0.4 Kg.
Power:	230V AC (5W).
Equipment Classification:	Class 2 Device, Continuous Operation.
Water protection:	Not protected.
Alarm:	Audio Visual.
Alarm range:	Default low = 18, High = 75. User selectable low = 15 TO 90. User selectable High = 21 TO 100.
Display:	LCD.